

Gnotobiotic Transfer Chamber

CRL Part # WPGTC2x2x3



At left – A typical Gnotobiotic Transfer Chamber with mechanically sealed doors on either end allowing the ability to transfer items from one environment to another eliminating cross contamination.

Note: The access doors can be supplied opaque or transparent with or without the access port shown.

Note: The access port shown comes with a polyurethane cap.

At right – Is our Gnotobiotic Transfer Chamber with dimensions as shown. However all of these dimensions can be customized to fit your application.

Dimensions:

Height: 24 inches
Length: 24 inches
Width: 32 inches



Gnotobiotic Transfer Chamber
CRL Part # WPGTC2x2x3